

PROTOCOL

CELL LYSIS FOR DOWNSTREAM ANALYSIS

RECOMMENDED MATERIALS AND REAGENTS

- Cells cultured with VitroGel system
- VitroGel Cell Recovery Solution (Cat# MS03-100)
- DPBS (no calcium, no magnesium)
- Cell lysis buffer
- Conical tubes (15 mL or 50 mL) or Microcentrifuge tube
- Serological pipettes or Micropipette; Low retention pipette tips
- Lab probe sonicator
- 37 °C water bath or dry bath
- Centrifuge

PROTOCOLS: (using 96 well-plate, 50 μ L gel/well as an example)

Protocol 1: Cell lysis with hydrogel together

1. Remove the medium covering the top of the hydrogel and wash the hydrogel three times with DPBS.
2. Add 100 μ L iced cold lysis buffer to hydrogel and mix the gel and buffer by pipetting up and down.
3. Transfer the mixture to a microcentrifuge tube and incubate on ice for 15 minutes.
4. Sonicate the mixture by using a lab probe sonicator: sonicate 2-5 seconds and rest for one minute on ice, repeat three times.
5. Incubate the mixture an additional 15 minutes (optional).
6. Centrifuge at 13,000 x g for 5 minutes at 4 °C.

Protocol 2: Harvest cells from hydrogel before cell lysis

1. Warm the VitroGel Cell Recovery Solution to 37 °C.
2. Take the cells out of the incubator and remove the medium covering the top of the hydrogel. Wash the hydrogel two times with DPBS.
3. Add 100 μ L warm VitroGel Cell Recovery Solution and mix the gel and solution by pipetting up and down 5-10 times.
4. Transfer the mixture to a microcentrifuge tube containing the 0.5-1 mL warm VitroGel Cell Recovery Solution.
5. Gently pipette the mixture up and down 5-10 times and incubate at 37 °C for 2-3 min. Repeat 3 times.
6. Centrifuge at 100 x g for 3-5 minutes at room temperature to collect the cell pellet (optimize the speed and time of centrifuge according to different cell types).

Optional: If there is still some hydrogel on top of the cell pellet, resuspend the cell with 1 mL warm cell recovery solution and repeat steps 5 and 6 one more time.



PROTOCOL

CELL LYSIS FOR DOWNSTREAM ANALYSIS

7. Remove the supernatant and add 100 μ L iced cold lysis buffer to the cell pellet.
8. Resuspend the cells in lysis buffer and incubate on ice for 15 minutes.
9. Sonicate the mixture by using a lab probe sonicator: sonicate 2-5 seconds and rest for one minutes on ice, repeat three times.
10. Incubate the mixture an additional 15 minutes (optional).
11. Centrifuge at 13,000 $\times g$ for 5 minutes at 4 $^{\circ}$ C.



Corporate Headquarters

TheWell Bioscience Inc.

675 US Highway 1

North Brunswick, NJ 08902

For research use only

Page 2 of 2

US TF: +1 866.3D.CELLS TEL: +1 973.855.4955 | support@thewellbio.com | www.thewellbio.com