

Ready-to-use, xeno-free hydrogel for 3D and 2D neuron culture



Xeno-free and biofunctional

100% synthetic. Animal & human origin-free hydrogel. No growth factors or proteins in the gel solution.



Supports 3D and 2D neuronal cell culture

Suitable for 3D and 2D culture of immortalized neuronal neuroblasts, iPSC-derived NSCs, neuron, and culturing of NSCs) as single cells or neurospheres.



Stable at room temperature

Easy-to-use system with room temperature operation. No ice bucket or chilled pipette tips required.




Injectable hydrogel for cell therapy

Its unique injectable properties make it an effective carrier for neuron-based cell therapy, spinal cord injury repair, and other regenerative medicine applications.



VitroGel® NEURON is the first-of-its-kind synthetic, xeno-free hydrogel specifically engineered to support both 2D and 3D neuronal culture and *in vivo* delivery. This groundbreaking platform eliminates the variability of animal-derived ECMs, such as Matrigel®, and provides a defined, reproducible, and scalable environment for neuronal growth, differentiation, and network formation.



2D Thin-Coating with VitroGel® NEURON Outperforms Animal-Based ECM



1:200

200X Dilution

4X more usage than animal-based ECM.




Long-Term Culture

Over 4 weeks of maintenance and post-differentiation culture.



Easy and Fast Operation

No need to remove the coating solution, just seed cells directly.



3D Neuron Culture Made Easy With VitroGel® NEURON



3D Neuron Differentiation

Supports both single cells and neuron spheroids.



Excellent Cell-Cell Communication

Enhanced cell extension and cell-matrix interactions.



Imaging-Friendly

Transparent and directly staining within the hydrogel for high-quality image data.

Data and References

2D Culture




Figure 1. VitroGel® NEURON 2D Thin Coating for Growth of iPSC-Derived NSCs and Confluence Analysis on VitroPrime™ Spread-Attach Plates. (A) iPSC-derived NSCs (CD34-eIPSC-NSC) were plated at 50,000 cells per well on 24-well VitroPrime™ Spread-Attach Plates coated with either VitroGel® NEURON (1:200) or Geltrex®, and cultured for 6 days. Images were captured using the Incucyte® S3 system. (B) NSC growth was quantified by percent confluence over time, with images collected every 8 hours and analyzed relative to the initial timepoint (T = 0 + 4 hrs).




Figure 2. 2D thin coating with VitroGel® NEURON hydrogel supports the development of Dopaminergic Neuron Progenitors (DNPs). (A) Images show robust survival of DNP colonies, and elaboration of neurites after 21 hours; (B & C). IncuCyte Neurotrack analysis software was used to quantitate the expansion of the area of DNP colonies, and their enhanced complexity of neurite elongation.

3D Culture




Figure 3. VitroGel® NEURON sustains long-term 3D neuronal differentiation. Bright field microscopy images depicting 3D B35 neuronal neuroblast differentiation in VitroGel® NEURON on days 7, 14, and 21. The images were obtained with the Keyence BZ-X microscope.




Figure 4. VitroGel® NEURON hydrogel supports the 3D development of dopaminergic neurons. iPSCs were 3D cultured in VitroGel® STEM with RocketCell™ IPSC Xeno-Free Growth Medium (Cat. No. RC02-GM) for 6 passages, and were differentiated into Floor Plate Neuro-Epithelial Spheroids. The spheroids were recovered from the hydrogel using VitroGel® Organoid Recovery Solution and embedded in VitroGel® NEURON.

(A) Brightfield images of spheroid on day 0 (A, left) that were differentiated for 16 days (A, right). (B) The samples were fixed and stained with beta-III-tubulin antibody and DAPI, depicting that the majority of the spheroids are neurons with the ability to form axonal connections. The cultures were visualized with the Keyence BZ-X microscope at 4X and 10X magnifications.

Learn more about
VitroGel® NEURON
thewellbio.com/VITROGEL-NEURON

Product	Cat No.	Size
VitroGel® NEURON	VHM07	10 mL
	VHM07S	2 mL

Complementary Products

RocketCell™ NSC Xeno-Free Growth Medium



Cell culture media engineered to promote NSC growth and maintenance.



Cat No: RC01-GM

VitroPrime™ Spread-Attach Plate




Supports exceptional 2D and 3D neuron growth. No PDL-coated plates required.



Plates Available: 6-well, 12-well, 24-well, 48-well, 96-well

CytoGrow™ Growth Factors



Premium quality growth factors with research-friendly price for advanced 3D/2D cell culture and beyond.



Scan the QR code to view the complete list of growth factors.