

## CytoGrow™ Recombinant Anti-Human CD52 mAb (Human)

Catalog Numbers: CG012-A, CG012-B, CG012-C

Revision 1.0

### Specifications

|                       |                                                                                                    |
|-----------------------|----------------------------------------------------------------------------------------------------|
| Synonym               |                                                                                                    |
| Species               | Human                                                                                              |
| Size                  | 50µg/1mg                                                                                           |
| Tag                   | Tag Free                                                                                           |
| Purity                | ≥95%, by SDS-PAGE (under reducing (R) & Non-reducing conditions, visualized by Coomassie staining) |
| Endotoxin             | ≤1 EU/mg by the LAL method                                                                         |
| Expression System     |                                                                                                    |
| Expression Region     | Mouse IgG2a/kappa                                                                                  |
| pH                    | 7.0-8.0                                                                                            |
| Appearance            | Colorless clear liquid                                                                             |
| Formulation           | 0.22 µm-filtered solution containing PBS, pH 7.4                                                   |
| Shipping              | The product is shipped with blue ice                                                               |
| Storage and Stability | Store in a manual defrost freezer and avoid repeated freeze-thaw cycles.                           |
|                       | Please keep under sterile conditions.                                                              |
|                       | 24 months at -20°C to -80°C<br>6 months at 2°C to 8°C                                              |