

CytoGrow™ Recombinant IL-23 (Human)

Catalog Numbers: CG50-A, CG050-B, CG050-C

Revision 1.0

Specifications

Synonym	IL-23 p19/IL-12 p40; IL23; IL-23A; Interleukin 23; SGRF
Species	Human
Size	10µg/50µg/1mg
Tag	Tag Free
Purity	≥95%
Endotoxin	≤10 EU/mg
Expression System	CHO
Expression Region	Ile23-Ser328 & Arg23-Pro189
pH	7.0-8.0
Appearance	
Formulation	Lyophilized
Shipping	The product is shipped with blue ice
Reconstitution	It is recommended to redissolve in sterile deionized water.
Storage and Stability	Lyophilized state: 36 months at -20°C to -80°C After reconstitution and under sterile conditions. 6 months at -20°C to -80°C 7-10 days at 2°C to 8°C Store in a manual defrost freezer and avoid repeated freeze-thaw cycles.

Reconstitution

1. Reconstitute protein to concentrations of 100 ug/mL or higher.
2. Use sterile deionized water to solubilize the lyophilized protein.
3. Gently tap glass vial to help collect the powder to the bottom.
4. Add desired amount of buffer, observe after 2 min. Gentle agitation by hand may be required.
5. Prepare aliquots of solution. Freeze at -20°C to -80°C to avoid repeated freeze thaw cycles.