

CytoGrow™ Recombinant LR3 IGF-1 (Bovine)

Catalog Numbers: CG055-A, CG055-B, CG055-C

Revision 1.0

Specifications

Synonym	
Species	Bovine
Size	10µg/50µg/1mg
Tag	A C-terminal 6-His tag
Purity	≥95%, by SDS-PAGE (under reducing (R)&Non-reducing conditions,visualized by Coomassie staining)
Endotoxin	≤10 EU/mg by the LAL method
Expression System	HEK293
Expression Region	G50-A119
pH	7.0-8.0
Appearance	White powder,Colorless clear liquid after reconstitution
Formulation	Lyophilized from a 0.2 µm filtered solution containing PBS,5% mannitol and 0.01% Tween 80, pH 7.4
Shipping	The product is shipped with blue ice
Reconstitution	It is recommended to redissolve in sterile deionized water.
Storage and Stability	Lyophilized state: 36 months at -20°C to -80°C After reconstitution and under sterile conditions. 6 months at -20°C to -80°C 7-10 days at 2°C to 8°C Store in a manual defrost freezer and avoid repeated freeze-thaw cycles.

Reconstitution

1. Reconstitute protein to concentrations of 100 ug/mL or higher.
2. Use sterile deionized water to solubilize the lyophilized protein.
3. Gently tap glass vial to help collect the powder to the bottom.
4. Add desired amount of buffer, observe after 2 min. Gentle agitation by hand may be required.
5. Prepare aliquots of solution. Freeze at -20°C to -80°C to avoid repeated freeze thaw cycles.